






C R E A T I V E L I V E

# Charles Krebs


Beyond Macro Photography: Into the Microscopic World























# Transition from Macro to Microphotography

## Beyond Macro Photography

- Macro lenses can be focused to .5X or 1X
- A variety of techniques can be used between 1X and 5X
- For magnifications of about 5X and higher, the best results are obtained utilizing microscope objectives



# Into the Microscopic World

- Difficulties encountered in close-up photography become extreme
- Optics, and some necessary equipment not familiar to most photographers
- Focus stacking becomes necessary in many cases
- The precision of a true microscope greatly facilitates working at very high magnifications

7





magnification = image size / object size

magnification = recorded width / sensor width

magnification = 4.46 / 22.3

magnification = 5X



## Between 1X and 5X

- Supplementary close-up lenses
- Stacked lenses
- Reverse mounted lenses
- Canon 65mm f/2.8 MP-E, 1X-5X lens












































2011

D



2011

D




2020















# The Three Main Optical Issues

- Depth-of-field
- Diffraction
- Vibration

# Depth-of-field

- As magnification increases, depth-of-field decreases very dramatically
- Smaller apertures can increase DOF, but increases diffraction and can significantly lower resolution
- Careful subject orientation will make best use of available DOF
- Combining the sharpest portions of images taken at various focus depths, or “focus stacking”, is a very effective way to create images of 3-dimensional subjects at high magnifications







Single image



86 image stack







# Moth Eggs



Single image




291 image stack























Single image



332 image stack



Mag (objective)	micron	mm	inch
1X (macro lens)	1000	1	0.039370
5X (5/0.15)	30	0.03	0.001181
5X (lens at f/2.8)	30	0.03	0.001181
10X (10/0.28)	8	0.008	0.000315
20X (20/0.40)	3.4	0.0034	0.000134
50X (50/0.50)	2.2	0.0022	0.000087
50X (50/0.80)	0.86	0.00086	0.000034

A standard CD is 1.2mm (1200 microns) thick.

A pin-head has a diameter of about 1.0 mm (1000 microns).

A typical grain of salt is approximately 0.3 mm (300 micron) on a side.

A red blood cell is 6-8 micron in diameter

A sheet of 20lb copy paper is 97-114 microns thick.



# Optical issues

- Depth-of-field
- Diffraction
- Vibration

# Diffraction

- Diffraction is caused by the wave nature of light, and in photography is directly related to **effective aperture** size
- As diffraction effects increase, the resolution and contrast of a photographic image is reduced
- Common f/#'s that we are familiar with are only valid when the lens is focused at “infinity” ( $\infty$ )
- The actual working, or *effective aperture*, becomes smaller when a lens is focused closer by extension




## Effective aperture ( $f_{\text{eff}}$ )


$$f_{\text{eff}} = f (m+1)$$

where  $f$  = f-number marked on lens,  $m$  = magnification on sensor  
(Note: This is a basic equation for symmetrical lenses, and is not valid for all lenses)

5X on APS-C sensor





17 image stack  
lens at  $f/3.5$   
effective aperture  $f/21$



Single image  
lens at  $f/16$   
effective aperture  $f/96$





# Effective aperture benchmark for sensor size

- Micro 4/3:  $f_{\text{eff}}/16-22$
- APS-C:  $f_{\text{eff}}/22$
- 35mm:  $f_{\text{eff}}/32$

At these benchmarks diffraction has already had some effect on image quality, but they reflect a reasonable compromise given available optical solutions, and average image usage. *Your image uses may warrant a larger or smaller value.*



# Optical Issues

- Depth-of-field
- Diffraction
- Vibration

# Vibration

- SLR mirror mechanism
- Mechanical shutter
- Handling. After a camera control has been touched, or the focus moved to a new point, equipment vibrations must be given some time to settle out
- External, or environmental, sources must be assessed



# Camera shutter vibration

- A fully mechanical shutter will negatively impact image quality
  - Use long shutter speeds from 1 to 4 seconds
  - Use electronic flash
    - Use long exposure with second curtain sync electronic flash
- Some cameras now have electronic first shutter curtain (EFSC)
- A few cameras now have fully electronic shutter option

# Choosing Microscope Optics

- Most conventional lenses used for macro have a maximum f-number of about  $f/2.8$ , but typically need to be closed down about 1-1.5 stops for the best quality
- As a result, if we want an effective aperture no smaller than  $f_{\text{eff}}/22$ , our useful magnification with that optic is in the range of 3X to 5X
- Microscope objectives can have large light gathering ability that provides a large effective aperture, so resolution is higher than with conventional macro optics



# Microscope objectives

- Light gathering ability denoted by **Numerical Aperture (NA)**, not f-number
- NA... key to resolution, analogous to maximum aperture
- Fixed magnifications such as 4/5X, 10X, 20X, 40/50X, 60X, 100X
- Designed to be used on a microscope, some can be used on bellows/tubes/lenses

## 10/0.28 microscope objective

Roughly equivalent to a 20mm f/1.6 lens

At 10X, Effective Aperture = f/18

## Regular lens at 10X

Marked f/2.8, actual Effective Aperture = f/31

Marked f/4, actual Effective Aperture = f/44

---

## 20/0.40 microscope objective

Roughly equivalent to a 10mm f/1.2 lens

At 20X, effective Aperture = f/24


## Regular lens at 20X

Marked f/2.8, actual Effective Aperture = f/59


Marked f/4, actual Effective Aperture = f/84



## Finite Type Objective



## Infinity Type Objective



# Important objective characteristics

## Degree of correction

- Different designed levels of optical and color correction... plan and non-plan, achromat, semi-apochromat, and apochromat

## Biological and Industrial objectives

- Biological are historically intended for transmitted lighting methods, and typically with a coverglass covered subject
- Industrial (or metallurgical) are primarily intended for opaque subjects that are not covered by a coverglass



# Working Distance (WD)

- The distance between the tip of the objective and the subject
- Highly variable and often far too small to be used off of a microscope with subjects using reflected light.
- WD typically not a major concern on a biological microscope

# Limited image circle

- Microscope objectives produce a relatively small image circle
- 18-20mm diameter image circle for older ones, 25-27mm diameter for newer super wide-field designs
- Pay attention to sensor size/image circle size. In many cases the actual usable image is significantly larger than specified



# Need for corrective eyepieces on some finite objectives

- Many older objectives used chromatically corrective eyepieces and photoeyepieces as the last stage of color correction. Avoid these for use with bellows/tubes
- No current infinity corrected optical systems use chromatically corrective eyepieces
- Two brands, Leica and Zeiss, currently utilize color corrective tube lenses. Best to avoid these mounted on a camera lens or bellows/tubes














# Microscope Optics Setups



Mitutoyo  
M26 x 0.7



Nikon CFI60  
M25 x 0.75




RMS thread  
M20.32





















## Mitutoyo M Plan Apo, LWD series, infinity type

5/0.15	WD = 34mm
7.5/0.21	WD = 35mm
10/0.28	WD = 33.5mm
20/0.46	WD = 20mm

## Olympus Plan Fluorite, infinity type

MPLFLN 5/0.15	WD = 20mm
MPLFLN 10/0.30	WD = 11mm
LMPLFLN 20/0.4	WD = 12mm
LMPLFLN 50/0.5	WD = 10.6mm

## Nikon CF M Plan Achromat 210/0 series, finite type

5/0.10 WD = 20mm (APS-C and smaller)

10/0.25 WD = 9mm

20/0.40 LWD WD = 6mm

40/0.50 ELWD WD = 10.1

60/0.70 ELWD WD = 4.9mm

## Nikon CFN series, finite 160mm

CFN Plan Apo 4/0.2 WD = 15mm (APS-C and smaller)

CFN Plan Achro 4/0.13 WD = 16.2mm (APS-C and smaller)

CFN Plan Achro 10/0.30 WD = 9.2mm (APS-C and smaller)



Noteworthy for low price, good performance, new objective

Nikon CFI BE Plan Achromat 4/0.01 (MRN70040, WD=25) infinity

Nikon CFI BE Plan Achromat 10/0.25 (MRN70100, WD=6.7mm) infinity

Nikon CFI60 Plan Achromat 10/0.25 (MRL00102, WD=10.1mm) infinity

# Focusing Platforms

- Must move accurately, repeatedly, in very fine increments
- Increments need to be less than DOF
- Minimal “wobble” around movement axis







Newport model 433 translation stage





































UMPlanFI N  
2.50  
26.5

UISA

LMPFLN50X

UMFLN10X


UISA

UMPlanFI  
5X/0.15  
JAPAN










# The Microscope

- Provides the ability to explore, study, and photograph the smallest of subjects
- Complete extremely high precision set-up solution with extremely fine movements for focus, x/y stage
- Can readily use high very NA, very small working distance, objectives
- Lighting techniques and methods that are very hard or impossible to create with a free standing, open set-up










































































































# Focus Stacking

- A series of images are taken, progressively shifting the focus point slightly with each image
- The focus shift should be less than the DOF (about 75% of DOF)

5/0.15	30 micron steps
10/0.28	8 micron steps
20/0.4	3.5 micron steps
50/0.5	1.5 micron steps



# Focus Stacking

- Camera should be set to a preset white balance (*not* AWB)
- Camera should be set for manual exposure
- Flash should be used at a manual setting
- Diffused lighting will generally give best results
- Be sure to use the techniques needed to minimize vibration problems

# Focus Stacking Software

- Helicon Focus
- Zerene Stacker
- CombineZP (free)
- Enfuse used with Hugin (free)



C R E A T I V E L I V E

# Charles Krebs

Beyond Macro Photography: Into the Microscopic World

Bonus Materials

To find the aperture size of a microscope in terms of a photographic “f-stop”, use the following relationship. This allows you to compare it to the maximum f-number of a conventional lens:

$$\text{f-number} = 1/(2 \cdot \text{NA}) * m/(m+1)$$

where m is magnification

To find the effective aperture of a microscope objective when used at its magnification you can use the following relationship:

$$\text{Effective aperture} = m/(2 \cdot \text{NA})$$



*Using the formulae on the previous page provides the following:*

### 10/0.28 microscope objective, at 10X

Maximum aperture (for comparison purposes) =  $f/1.6$

Effective Aperture =  $f/18$

### Regular lens such as macro lens, enlarging lens:

Marked  $f/2.8$ , at a 10X magnification,  $f_{\text{eff}} = f (m+1)$

Effective Aperture =  $f/31$

Marked  $f/4$ , at a 10X magnification,  $f_{\text{eff}} = f (m+1)$

Effective Aperture =  $f/44$

The microscope objective can exhibit much higher resolution than the conventional lens because the losses from diffraction are much less

*At 20X, the resolution advantage of a microscope objective is extremely significant*

### 20/0.40 microscope objective, at 20X

Maximum aperture (for comparison purposes) =  $f/1.2$

Effective Aperture =  $f/24$

### Regular lens such as macro lens, enlarging lens:

Marked  $f/2.8$ , at a 20X magnification,  $f_{\text{eff}} = f(m+1)$

Effective Aperture =  $f/59$

Marked  $f/4$ , at a 10X magnification,  $f_{\text{eff}} = f(m+1)$

Effective Aperture =  $f/84$

The microscope objective can exhibit much higher resolution than the conventional lens because the losses from diffraction are significantly less



A superb resource for photographers wishing to learn more about this field of photography is:

[www.photomacrography.net](http://www.photomacrography.net)